

In ImageJ the equivalent to ‘Redo’ is the **Process**▷Repeat Command [R], that re-runs the previous used command (skipping **Edit**▷Undo [z] and **File**▷Open... [o] commands).

SEE ALSO: **Plugins**▷Utilities▷Reset... , [Multi Undo](#) plugin

## 7 Image Types and Formats

Digital Images are two-dimensional grids of pixel intensities values with the width and height of the image being defined by the number of pixels in  $x$  (rows) and  $y$  (columns) direction. Thus, pixels (picture elements) are the smallest single components of images, holding numeric values – pixel intensities – that range between black and white. The characteristics of this range, i.e., the number of unique intensity (brightness) values that can exist in the image is defined as the bit-depth of the image and specifies the level of precision in which intensities are coded, e.g.: A 2-bit image has  $2^2 = 4$  tones: 00 (black), 01 (gray), 10 (gray), and 11 (white). A 4-bit image has  $2^4 = 16$  tones ranging from 0000 (0) to 1111 (16), etc. In terms of bits per pixel (bpp), the most frequent types of images (**Image**▷Type▷) that ImageJ deals with are ([ImageJ2 supports many more types of image data](#)):

- 8-bit** Images that can display 256 ( $2^8$ ) gray levels (integers only).
- 16-bit** Images that can display 65,536 ( $2^{16}$ ) gray levels (integers only).
- 32-bit** Images that can display 4,294,967,296 ( $2^{32}$ ) gray levels (real numbers). In 32-bit images, pixels are described by [floating point](#) values and can have ANY intensity value including *NaN* (Not a Number).
- RGB Color** Color Images that can display 256 values in the Red, Green and Blue channel. These are 24-bit ( $2^{3 \times 8}$ ) images. RGB color images can also be 32-bit color images (24-bit color images with additional eight bits coding alpha blending values, i.e., transparency).

### Native Formats

Natively (i.e. without the need of third-party plugins) ImageJ opens the following formats: **TIFF**, **GIF**, **JPEG**, **PNG**, **DICOM**, **BMP**, **PGM** and **FITS**. Many more formats are supported with the aid of plugins. These are discussed in Non-native Formats.

- TIFF** (Tagged Image File Format) is the ‘default’ format of ImageJ (cf. **File**▷Save [s]). Images can be 1-bit, 8-bit, 16-bit (unsigned<sup>1</sup>), 32-bit (real) or RGB color. TIFF files with multiple images of the same type and size open as Stacks or Hyperstacks. ImageJ opens lossless compressed TIFF files (*see* II Image Types: Lossy Compression and Metadata) by the LZW, PackBits and ZIP (Deflate/Inflate) [2] compression schemes. In addition, TIFF files can be opened and saved as ZIP archives. Tiff tags and information needed to import the file (number of images, offset to first images, gap between images) are printed to the Log Window when ImageJ is running in *Debug Mode* (**Edit**▷Options▷Misc... , *see* Settings and Preferences).

---

<sup>1</sup>A numeric variable is signed if it can represent both positive and negative numbers, and unsigned if it can only represent positive numbers.

- DICOM** (Digital Imaging and Communications in Medicine) is a standard popular in the medical imaging community. Support in ImageJ is limited to uncompressed DICOM files. DICOM files containing multiple images open as Stacks. Use **Image**▷**Show Info...** [i] to display the DICOM header information. A DICOM sequence can be opened using **File**▷**Import**▷**Image Sequence...** or by dragging and dropping the folder on the 'ImageJ' window. Imported sequences are sorted by image number instead of filename and the tags are preserved when DICOM images are saved in TIFF format. ImageJ supports custom DICOM dictionaries, such as the one at [http://imagej.nih.gov/ij/download/docs/DICOM\\_Dictionary.txt](http://imagej.nih.gov/ij/download/docs/DICOM_Dictionary.txt). More information can be found at the [Center for Advanced Brain Imaging](#).
- FITS** (Flexible Image Transport System) image is the format adopted by the astronomical community for data interchange and archival storage. Use **Image**▷**Show Info...** [i] to display the FITS header. More information [here](#).
- PGM** (Portable GrayMap), **PBM** (Portable BitMap) and **PPM** (Portable PixMap) are simple image formats that use an ASCII header. More information [here](#).
- AVI** (Audio Video Interleave) is a container format which can contain data encoded in many different ways. ImageJ only supports uncompressed AVIs, various YUV 4:2:2 compressed formats, and PNG or JPEG-encoded individual frames. Note that most MJPG (motion-JPEG) formats are not read correctly. Attempts to open AVIs in other formats will fail.

SEE ALSO: Non-native Formats, II Image Types: Lossy Compression and Metadata, X Warning on JPEG Compression

## Non-native Formats

When opening a file, ImageJ first checks whether it can natively handle the format. If ImageJ does not recognize the type of file it calls for the appropriate reader plugin using [HandleExtraFileTypes](#), a plugin bundled with ImageJ. If that fails, it tries to open the file using the [OME Bio-Formats library](#) (if present), a remarkable plugin that supports more than [one hundred of the most common](#) file formats used in microscopy. If nevertheless the file cannot be opened, an error message is displayed.

Because both these plugins are under active development, it is important that you keep them updated. The OME Bio-Formats library can be updated using its self-updating plugin (**Plugins**▷**LOCI**▷**Update LOCI Plugin...**) or Fiji's built-in updater (**Help**▷**Update Fiji...**). The following websites provide more information on the OME Bio-Formats:

- <http://loci.wisc.edu/bio-formats/imagej>
- <http://fiji.sc/Bio-Formats>
- <http://loci.wisc.edu/bio-formats/using-bio-formats>

In addition, the ImageJ web site lists [more than sixty plugins](#) that recognize more 'exotic' file formats. The ImageJ Documentation Portal also maintains a (somewhat outdated) [list of file formats](#) that are supported by ImageJ.

SEE ALSO: Native Formats, **File**▷**Import**▷, II Image Types: Lossy Compression and Metadata, X Warning on JPEG Compression, [Acquisition plugins](#), [Input/Output plugins](#)